


Inlet Supplies for PerkinElmer GCs With CAP Inlet



Septa (11 mm) for PerkinElmer GCs

| Thermolite Septa | |
|------------------------------|-------|
| Usable to 340 °C inlet temp. | |
| qty. | cat.# |
| 50-pk. | 27142 |
| 100-pk. | 27143 |

| Premium Non-Stick BTO Septa | |
|------------------------------|-------|
| Usable to 400 °C inlet temp. | |
| qty. | cat.# |
| 50-pk. | 27090 |
| 100-pk. | 27091 |

Silicone O-Rings

| Description | Max. temp. | Similar to PE part # | qty. | cat.# |
|---|------------|----------------------|--------|-------|
| Silicone O-Rings for PerkinElmer Auto SYS XL or Clarus w/CAP Injector | 250 °C | N6101374 | 10-pk. | 20262 |

Sky® Inlet Liner

| ID x OD x L | qty. | cat.# |
|--|-------|---------|
| Precision, Sky Technology, Borosilicate Glass with Quartz Wool | ea. | 23450.1 |
| 4.0 mm x 6.2 mm x 92.1 mm | 5-pk. | 23450.5 |

Capillary Ferrules for 1/16-Inch Compression-Type Fittings

| Ferrule ID | Fits Column ID | qty. | Vespe/Graphite (60/40) |
|------------|-----------------------------|--------|------------------------|
| 0.4 mm | 0.10/0.15/0.18/0.25/0.28 mm | 10-pk. | 20211 |
| 0.5 mm | 0.32 mm | 10-pk. | 20212 |
| 0.8 mm | 0.45/0.53 mm | 10-pk. | 20213 |



Restek GmbH

Schaberweg 23 * 61348 Bad Homburg

Tel: 06172 / 27 97 0 * Fax: 06172 / 27 97 77

info@restekgmbh.de * www.restekgmbh.de